Content of manuals:

EN - English

ES - Spanish

FR - French

DE - German

NL - Dutch

DA - Danish

SE - Swedish

NO - Norwegian

FI - Finnish

TR - Turkish

PO - Portuguese

IT - Italian

PainSensor

User manual

VERSION 16

ENGLISH

DHF-00069

PainSensor REF 1002

Part number 400102

Manufacturer/Distributor:

Med-Storm Innovation AS
Gimle Terrasse 4
NO-0264 Oslo
Norway

Telephone:

+47 90 93 98 10

Internet:

<http://www.med-storm.com/>

Trademark:

MedStorm Innovation AS

Date of publication:

12.01.2025

IMPORTANT

The user manual covers the operation of the MED-STORM PainSensor REF1002.

Read all instructions, warnings and precautions prior to use.

Only a trained physician or nurse should use the system.

Users of the equipment must be familiar with the medical aspects of the conditions for which the MED-STORM PainSensor is used. All users must complete the training program before using the PainSensor.

MED-STORM considers itself responsible for any effects on safety, reliability and performance of the equipment only if:

- assembly operations, extensions, re-adjustments, modifications or repairs are carried out by persons authorized by MED-STORM, and
- the electrical installation complies with national standards, and
- installation and configuration of software is carried out by persons authorized by MED-STORM, and
- no other software is installed on the computer or the PainSensor unless explicitly accepted by MED-STORM, and
- the equipment is used in accordance with the product documentation.

MED-STORM makes no warranty of any kind with regard to this material, including, but not limited to, the implied warranties of merchantability and fitness for a particular purpose.

MED-STORM shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material.

DISCLAIMER

THE MED-STORM PAIN SENSOR IS NOT A SUBSTITUTE FOR YOUR PROFESSIONAL JUDGMENT. MED-STORM SHALL NOT BE LIABLE IN ANY MANNER WHATSOEVER FOR THE RESULTS OBTAINED THROUGH THE USE OF THE PAIN SENSOR. PERSONS USING THE PAIN SENSOR ARE RESPONSIBLE FOR THE SUPERVISION, MANAGEMENT AND CONTROL OF THE PAIN SENSOR.

TABLE OF CONTENTS

TABLE OF CONTENTS.....	3
1. Introduction.....	4
1.1 Intended use.....	4
1.2 Normal use.....	4
1.3 Intended user.....	4
1.4 Indications for use.....	4
1.5 Contraindications for use.....	4
1.6 Pre use checks.....	4
2. Warnings, limitations, and residual risks.....	4
2.1 Electrical shock hazard.....	6
2.2 Warning and information symbols.....	6
3. System overview.....	7
3.1 PainSensor unit.....	8
3.2 Electrodes.....	10
4. Operating instructions.....	11
4.1 Data security and network.....	11
4.2 Installation and setup.....	12
4.3 Button and LEDs.....	13
4.4 Error conditions.....	14
4.5 Getting started.....	15
5. Care and maintenance.....	18
5.1 Lifetime.....	18
5.2 Support information.....	18
5.3 Cleaning.....	18
5.4 Scrapping instructions.....	19
6. Appendix A - Set-up of the PSM connectivity box to connect with the PainSensor.....	20
7. Appendix B - Connecting the PainSensor to a stand-alone PC with Bluetooth connection.....	20
7.1 The PainSensor Application.....	20
7.2 Connect to PainSensor.....	21
7.3 Description of the interface.....	22
7.4 Connect the tablet with a Philips patient monitor.....	25
7.5 Store and analyze data.....	26
8. Appendix C - Environmental and handling conditions.....	29
9. Appendix D - Technical specifications.....	30
10. Appendix E - Regulations and Directives.....	30
11. Appendix F - Electromagnetic compatibility.....	31
12. Appendix G - Pre use checklist.....	35

1. Introduction

This user manual corresponds to PainSensor REF1002

1.1 Intended use

The PainSensor is intended to determine a patient's sensitivity to pain, awakening from pain during anaesthesia, abstinence symptoms and nerve block.

1.2 Normal use

Normal use is 24 hours/day, 200 days a year. The electrodes do not allow continuous use for more than 48 hours, after 48 hours the electrodes (consumable components) have to be changed. The PainSensor can be reused after adequate checks and cleaning (see section 5).

1.3 Intended user

Only trained physicians or nurses shall use the system.

1.4 Indications for use

Indications for use are:

- patients undergoing anaesthesia;
- postoperative patients;
- patients in the intensive care units;
- premature infants;
- patients exposed to regional nerve block;
- patients with withdrawal symptoms / abstinence.

1.5 Contraindications for use

- The device shall not be used on patients with skin conditions which may affect skin conductance e.g. injury of the skin.
- More than one device shall not be used on patients with electrically sensitive life support systems (e.g. implantable pacemaker or defibrillator).
- The device shall not be used when the patient has an injury affecting the sympathetic skin nerves.

1.6 Pre use checks

Before using the device, the user must fulfil training either provided by MedStorm Innovation AS, or by someone who has received training from MedStorm Innovation AS.

Users are invited to use the Pre use checklist provided in appendix G of this manual.

2. Warnings, limitations, and residual risks

- WARNING (1):** Read the entire user manual before operating this PainSensor
- WARNING (2):** It is the responsibility of the user to ensure that any applicable regulations regarding the operation of the PainSensor are observed.
- WARNING (3):** The PainSensor should not be used as the only source of information for pain.
- WARNING (4):** Modification of the PainSensor is not permitted.
- WARNING (5):** The maximum temperature is 43°C, contact time more than 15.
- WARNING (6):** Do not use more than ONE PainSensor at the same time on patients with an implanted pacemaker or defibrillator. This may risk cardiac problems.
- WARNING (7):** The device shall not be used on patients with skin conditions which may affect skin conductance, e.g. injury of the skin below the electrodes or when

the patient has an injury affecting the sympathetic skin nerves. Moreover, local nerve blocks at the measuring area will affect the method.

WARNING (8): Neuromuscular blockers may block the index in preterm infants

Charging and battery

WARNING (9): The PainSensor cannot and must not be used on the patient during charging. During charging, the measurement is prevented by hardware. The power supply cable also has to be removed from the PainSensor while the PainSensor is used on the patient. Charging should take place at a different place than where the patient is located.

WARNING (10): Only use the PainSensor with the power supply and power cable delivered and specified by MedStorm Innovation AS.

WARNING (11): Position the device in such a way that you can disconnect the power plug from the device at any time when the device is charging.

WARNING (12): The PainSensor can be used whenever needed. However, the battery capacity limits the continuous use between charging, the battery capacity is at least 15 hours.

WARNING (13): Do not touch the power supply unit /USB connector and the patient at the same time

Use with accessories

WARNING (14): The PainSensor must be used together with dedicated accessories.

WARNING (15): Only use Bluetooth devices that have the required software installed to connect to the PainSensor.

WARNING (16): If information from the PainSensor is absent on one tablet or PainSensor, you can use a different tablet or PainSensor from MedStorm Innovation to show the information

WARNING (17): The same electrodes do not allow continuous use for more than 48 hours.

Environmental conditions

WARNING (18): Do not use, transport or store above or below the recommended environmental intervals in Appendix D.

WARNING (19): Do not immerse the PainSensor or cables in any liquid or allow liquid to enter plugs or connections. Do not use cables if connectors become wet.

WARNING (20): Avoid any spillage of liquids on the PainSensor.

WARNING (21): Simultaneous use of electro surgery may disturb the measurements made with the PainSensor for a period of time.

IT

WARNING (22) If the PainSensor is used with other devices in an IT network as an ME system, the commissioning person must verify that the system complies with the applicable safety requirements and regulations

The PainSensor has no "Essential performance". The output of the device is one of several indicators that may be used to determine patient status, thus the absence of performance does not lead to an unacceptable risk.

REPORTING:




Any serious incident that has occurred in relation to the device should be reported to Med-Storm Innovation AS and the competent authority of the Member State in which the user and/or patient is established.

2.1 Electrical shock hazard

There are exposed voltages inside the PainSensor. There are no user-serviceable parts inside. Do not open the PainSensor. Send to qualified personnel approved by Med-Storm Innovation for servicing.

2.2 Warning and information symbols

Table 2-1: Warning and safety information symbols used on the PainSensor

	Caution is necessary when operating the PainSensor.		
	Read the entire instructions for use before operating this PainSensor		Class II electrical equipment for protection against electric shock (safety connection to electrical earth not required).
	Dispose of this unit as normal electronic household waste per local regulations according to the WEEE Directive		TYPE BF APPLIED PART
	Serial number of device		TYPE B APPLIED PART
	Product model identifier in catalogue		The box (and contents) should be kept dry
	Batch designation		Date of manufacture
	Medical device		Manufacturer
	Unique device Identifier		Temperature limit
	Atmospheric pressure limitation		Humidity limitation

CE mark and notified body number

QR code for accessing the user-manual online

3. System overview

The PainSensor is an electronic conductance meter for detecting skin conductance changes on palmar and plantar skin sites to determine a patient’s sensitivity to pain. A sketch of the system is shown in Figure 3-1.

Note.

The use of accessories, transducers and cables other than those specified may result in increased emission or decreased immunity of the system.


Figure 3-1 System overview. The PainSensor with wristband, charger and electrodes are shown on the left side of the bluetooth symbol. The connectivity options are shown below the bluetooth symbol (tablet with PSS software) and on the right side of the bluetooth symbol (connectivity box to patient monitors).

The following parts of PainSensor are included in the package from MedStorm Innovation AS

Measurement equipment
PainSensor (8,2 x 2,6 x 5,4 cm)
Power supply power cable
Power supply unit for PainSensor (1.8 x 4,2 x 6 cm)
Power supply with internal adapter plugs
Wristband (2 x 0.1 x 15 cm)
Manual, English, Norwegian, German, Dutch, Swedish, Danish, Spanish, French
Firmware
Microcontroller

BLE

The following parts can be used together with the PainSensor:

Accessories
Electrodes for adult and children
Electrodes for infants
Electrodes for preterm infants
PSM connectivity box
External PC screen with power supply, power cable, and PainSensor Software (PSS) Application

Applied parts:

Description	Inside/Outside	High touch probability	Low touch probability	Applied part (Type)
Housing (+intermediate ring, bend relief and wristband)	Outside	x		Type B
Sticker (HMI)	Outside	x		Type B
Electrode cable	Outside	x		Type B
Power supply	Outside		x	
USB cable	Outside		x	
Electrodes	Outside	x		Type BF

Bluetooth transmitter specification:

Feature	Specification
Bluetooth®	V5.0 – Single mode, Concurrent master and slave, Diffie-Hellman based pairing
Frequency	2.402 - 2.480 GHz
Maximum Transmit Power Setting	+4 dBm
Minimum Transmit Power Setting	-20 dBm

3.1 PainSensor unit

The PainSensor is in the product category Pain Monitoring System. The PainSensor, Figure 3-2, has two connectors, one button, and two LEDs. The connector to attach the electrodes is one connector [1]. The power supply inlet [2] is the other connector, and both are placed on one side of the PainSensor. The LEDs [3][4] and the button [5] are found at the front of the PainSensor. One LED [3] shows the status of the PainSensor, and one LED [4] shows the charging status of the PainSensor. There is a plug covering the charging USB port when the charging cable is not attached.


Figure 3-2: PainSensor unit.

3.1.1 Power supply and power cable

Figure 3-3: PainSensor power supply

The power supply unit and power cable used with the PainSensor is medical grade and provides 2 MOPP.

Specifications:

Model	HDP12-MD05024U Certified to IEC 62368-1 and IEC 60601 standards. For Medical & ITE Applications.
Input Voltage	100 ~ 264 VAC
Frequency	50/60 Hz
Input Current	0.35A
Output Voltage	5.0V
Output Current (max)	2.4 A

Do not use any other than the recommended power supply unit and power cable (HDP12-MD05024U). The power supply unit and cable are shown in, Figure 3-3. The power supply unit has one input for

connecting the cable. The orange rectangles mark which side of the cable should be connected with the power supply input. To disconnect the PainSensor from the mains supply, the mains plug must be pulled out of the socket. Charging of the device should take place a different place than where the patient is located.

3.2 Electrodes

The system can use 3 types of electrodes supplied from Med-Storm Innovation AS. #6101 electrodes are for adults and children (Figure 34), whereas #6102 electrodes are for infants (Figure 35). #6103 electrodes are for preterm infants less than 37 weeks of gestational age (Figure 36). All infants shall have the bear under the sole of the foot. Ensure reading IFU for electrodes before use.

Figure 3-4: Adult electrodes

Figure 3-5: Infant electrodes

Figure 3-6: Preterm Infant electrodes

4. Operating instructions

4.1 Data security and network

The PainSensor is integrated into an IT network to transmit the measured and calculated data via Bluetooth Low Energy to a unit for display (as described in Appendix A and Appendix B). The network consists of a two-node network between the PainSensor and the receiving unit. Only Bluetooth devices that have the required software to read out the data can be used

The information flow is purely from PainSensor to the respective recipient.

No hazards arise if an IT network cannot provide the required features. A connection with the PainSensor is then not possible.

WARNING:

If the PainSensor is used with other devices in an IT network as an ME system, the commissioning person must verify that the system complies with the applicable safety requirements and regulations.

Attention:

The integration of the PainSensor into existing IT networks can lead to risks for the patient operator or for third parties which cannot be foreseen by the manufacturer.

The person commissioning the system should determine, analyse, evaluate and manage these risks to ensure safety.

The following changes to the IT network could lead to new risks and should therefore be analysed:

- Changes to the IT NETWORK configuration.
- Connecting additional elements to the IT NETWORK
- Removing elements from the IT NETWORK;
- "Update" of devices connected to the IT NETWORK.
- "Upgrade" of devices connected to the IT NETWORK.

Intended position of the operator is less than 0,4 m from the device. The following instructions describe all necessary steps required to set-up and operate the PainSensor.

Note.

Portable and mobile wireless communications equipment can affect medical electrical equipment.

Note.

Medical electrical equipment needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in Appendix F.

4.2 Installation and setup

Power supply:

1. Connect the power cable to the power supply unit (Only products delivered by MedStorm Innovation AS, check label).
2. Remove the plug covering the USB port for charging on the PainSensor, and connect the power cable to the PainSensor (Only products delivered by MedStorm Innovation AS, check label).
3. Connect the power supply unit to a wall socket

Connection of PainSensor to patient:

4. Make sure that the charger is disconnected from the PainSensor. The device has a battery and can, when fully charged, be used without connection to the power supply. The battery life enables at least for 15 hours of continuous measurement.
5. Place the electrodes on the patient and attach the electrodes to the device. *You will hear a click when it is well connected.* See description in section 4.2.1.
6. Attach the wristband with the PainSensor to the wrist/ankle of the patient, or place the PainSensor next to the patient.
 - 7a. For infants in the NICU, the PainSensor can be placed inside or outside the incubator, or on the wrist of the parents during skin-to-skin care. The sensor should not be strapped directly on the infant.
 - 7b. For use in the operating room, the PainSensor can be placed on the underside of the arm if that is considered more convenient.

4.2.1 Skin electrode placement

The electrodes can be attached to the patient with reliable measurement result for a maximum time of 48 hours.

4.2.1.1 Skin electrode placement on adults and children

The intended placement of the electrodes on adults and children is on the palm of the hand. The sole of the foot can be used in active patients or if the child is young. If using the sole of the foot, make sure that you place the electrodes in an area where the skin is not thick so the sweat reach the upper part of the skin. The skin does not need to be prepared before placing the electrodes.

Place the electrodes inside the palm of the hand or under the foot of the patient, see example in Figure 4-1. Place the electrodes such that the electrode in the middle is placed at the hypothenar eminence because this area on the palm gives highest stability and thus less movement artefacts.

Figure 4-1: electrode placement under the foot or in the palm of the hand

4.2.1.2 Skin electrode placement on infants

The intended placement of the electrodes on infants is under the foot. For infant electrodes (6102), the distance in between each electrode shall be at least 7 mm, with the middle electrode under the sole and the two others on each side of the angle. The electrode in the middle should be placed under the foot, as shown in Figure 4-2.

Figure 4-2: electrode placement on infants

4.2.1.3 Artefacts

Artefacts can be seen when moving the hand/foot where the electrodes are attached or by pulling at an electrode. If the measuring electrode is wrapped, e.g. with a bandage, the movement artefacts should be reduced/eliminated. It is therefore advised to wrap the electrodes if the patient is excessively moving the hand/foot where the electrodes are placed.

If the measuring electrode (middle electrode) is fastened to the extremity with regional block, no response to pain/noxious stimuli will be observed because the skin sympathetic nerves are blocked.

Artefacts can also be seen if the electrodes are attached to injured skin.

Artefacts have been seen in the registration during tetanic stimuli if there were more than one EEG (Electroencephalography) monitor connected to the patient.

4.3 Button and LEDs

The LEDs show a continuous light in different colours, with different meanings. These are depicted in Figure 4-3.

Continuous light	
 PainSensor is ready for use.	 Bluetooth connection successful
 Measure and transmit	 Charging
 Low battery	 Charging complete

Figure 4-3: Meaning of the different colors of the LED.

The button can be pressed long or short, as depicted in Figure 4-4.

- Short press (1 second or shorter)
 - Switch on device
- Long press (2 seconds or longer)
 - Switch off device

Figure 4-4: Functioning of the button press.


4.4 Error conditions

The PainSensor measures very small changes in skin conductance and is extremely sensitive. Simultaneous use of electro surgery will for example disturb the measurements made with the PainSensor.

If the PainSensor loses contact with the M electrode (middle), if the microsiemens level is below 1 for a certain time, or if the patient is moved or moving creating movement artifacts (the PainSensor has to be worn by the patient to detect movement artifacts), the system will automatically recognize this and indicate this with a “bad signal quality” error on the PC/monitor. If using the software program on a tablet, the indices will be replaced by “- -” instead of a number and the trend graph will show vertical lines to indicate “bad signal” in the PainSensor Application. When using the PSM Connectivity Box, bad signal will be indicated with a “1” while if the bad signal is not present this will be indicated with a “0”. If the patient has severe eudema, the fluid in the skin may induce pressure on the sweat glands, and the sweat will not reach the surface of the skin.

If “bad signal” is indicated the following should be checked:

- a) Are the electrodes properly connected to the PainSensor?
- b) Are the electrodes properly attached to the patient?
- c) Is the patient and/or PainSensor being moved?
- d) Is the micro Siemens level below 1?
- e) Does the patient have eudema?

*note that the software might have a slightly different design compared to this picture for different versions, but the functionality is the same.

4.5 Getting started

Note.

Ensure that the PSM connection box / stand-alone PC are prepared for connecting to the patient monitors. All information about the installing and operating the PSM connection box is found in Appendix A. For transmitting data to a stand-alone PC, please see Appendix B.

- a) Make sure that the PainSensor is charged.
- b) After the PainSensor is connected to the patient by following steps 4-7 in section 4.2, switch on the device. This is done by pressing the button one time. When the device is ready for use, a green light is shown on the PainSensor.
- c) In order to visualize the data, press “connect sensor” in the PSS software program on your tablet, or the “on” button on the PSM connectivity box. See appendixes for detailed information about this connectivity. The PainSensor will show a blue light when it is connecting to a device for visualizing the data.
- d) The PainSensor will then automatically connect to the chosen connectivity option, and will automatically start measuring and transmitting data, seen with the cyan light on the PainSensor.

It is possible to see that the device is powered either by detecting a green LED, a blue LED showing that the device successfully connected to Bluetooth, or when the cyan LED is constant for measuring.

- In case of failure to connect to bluetooth, the device will automatically reattempt to connect. If connection is not successful after some time, the PainSensor will switch off. Reattempt to connect by starting the process again from step 2 in this section.
- If you are unable to connect the device, try to use a different tablet / PSM connectivity box. If you are still unable to establish a connection, turn off the device by holding the button for a few seconds, remove the device from patient, and contact the manufacturer MedStorm Innovation at support@med-storm.com.

- e) Switch off the device by a long press of the button if you are done measuring, or if an error occurs.

4.5.1 Operating theatres and Intensive care units (sedated patients)

In operating theatres, “peaks/second” and “area under curve” are used. If the “Peaks per sec” measure is 0 the patient is sufficiently or too sedated. If there are peaks, 2 or more in the analysing window of 15 sec, there are bursts in the sympathetic nerves and the patient’s sensitivity to pain is reached during anaesthesia. The PainIndex is shown in Figure 4-5. Once the patient has emerged from general sedation, it is recommended to switch to "post operative" setting in the application.

Pain Index	Peaks per second Mirroring how often the skin sympathetic nerves are firing	Indication - – During Surgical stimuli / Procedural painful stimuli
0	0.00 - 0.06	No more analgesia is needed
1	0.07 - 0.12	Analgetic light and adequate level – the goal
3	0.13 - 0.19	More analgesia may be needed

5	0.20 - 0.26	More analgesia is needed
7	0.27 - 0.32	More analgesia is needed
8	0.33 - 0.39	More analgesia is needed
10	0.40 or higher	More analgesia is needed

If **Figure 4-5: The PainIndex, the corresponding peaks/second for each value, and the indication** the “area under the curve” increases to 100 the patient is about to wake up from the stimulus, and the patients possibly need more analgesics and hypnotics (forceful bursts in the sympathetic nerves when the patient is about to wake up). *Awakening index is a response that appears only if there is a strong nociceptive stimulus, it will not be triggered in case there is too low level of hypnotics and the anti-nociceptive drug level is adequate.* The awakening index is shown in Figure 4-6.

Awakening Index	AUC (unit: siemens-seconds) Mirroring how forceful the skin sympathetic nerves are firing	Indication - – During Surgical stimuli / Procedural painful stimuli
0	0.00 - 1.99	No more analgesia needed
40	2.00 - 4.99	More analgesia is possibly needed
100	5.00 – 10	More analgesia is needed

Figure 4-6: The awakening index, the corresponding area under the curve, and the indication. The awakening index is so far only validated in adults.

4.5.2 Post-operative and ICU (awake patients)

The measurement values shown in “post operative” mode are shown in Figure 4-7.

Pain Index	Peaks per second Mirroring how often the skin sympathetic nerves are firing	Indication - during ongoing painful stimuli
0	0.00 - 0.06	No more analgesia is needed
1	0.07 - 0.12	No more analgesia is needed
3	0.13 - 0.19	No more analgesia is needed
5	0.20 - 0.26	More analgesia may be needed
7	0.27 - 0.32	More analgesia is possible needed
8	0.33 - 0.39	More analgesia is probably needed
10	0.40 or higher	More analgesia is needed

Figure 4-7: Pain Index for post-operative mode and ICU mode for awake patients, the corresponding value, and the indication

4.5.3 Infants

In the “Infant” application mode the measurement values shown are shown in Figure 4-8.

Pain Index	Peaks per second Mirroring how often the skin sympathetic nerves are firing	Indication - during painful events
0	0.00 - 0.06	No more analgesia is needed
1	0.07 - 0.13	No more analgesia is needed

3	0.14 - 0.19	More analgesia may be needed
5	0.20 - 0.26	More analgesia is possibly needed
7	0.27 - 0.32	More analgesia is probably needed
8	0.33 - 0.39	More analgesia is needed
10	0.40 or higher	More analgesia is needed

Figure 4-8: The pain values for infants, the corresponding values and the indication

Note: The influence from muscle relaxant agents and atropine has not been validated for infants.

4.5.4 During abstinence symptoms

The PainSensor index when tailoring the need of analgesia during abstinence symptoms is shown in Figure 4-9.

Pain Index	Peaks per second <small>Mirroring how often the skin sympathetic nerves are firing</small>	Indication
0	0.00 - 0.06	No more analgesia is needed
1	0.07 - 0.13	No more analgesia is needed
3	0.14 - 0.19	More analgesia may be needed
5	0.20 - 0.26	More analgesia is possibly needed
7	0.27 - 0.32	More analgesia is probably needed
8	0.33 - 0.39	More analgesia is needed
10	0.40 or higher	More analgesia is needed

Figure 4-9: Pain index during abstinence symptoms. The color is visible when a PC is used to show the index. On the monitors, the values on the right hand side are shown

The abstinence index is so far only validated in infants.

4.5.5 Nerveblock

The PainSensor indices to assess the effect of regional peripheral nerve block are shown in Figure 4-10.

NerveBlock Index	Peaks per second <small>Mirroring how often the skin sympathetic nerves are firing</small>	Indication
0	0.00 - 0.06	0 for 2 minutes, the nerve is probably blocked
1	0.07 - 0.13	
3	0.14 - 0.19	
5	0.20 - 0.26	
7	0.27 - 0.32	
8	0.33 - 0.39	
10	0.40 or higher	

Figure 4-10: The NerveBlock index to assess the effect of regional peripheral nerve block.

The nerve block is possible to use in awake patients. If the patient is sedated, the nerves will already be blocked on spinal level. Secure that the index changes for a value higher than 0 peaks per sec before the nerve block to a value 0 per sec when the nerve block starts to work. The value should stay on 0 for at least 2 minutes to secure that nerve block has been obtained.

5. Care and maintenance

Routinely inspect all electrical plugs and connections of the medical device. Check that none of the cables are destroyed, check that the LED is working, check that the button is functioning. Do not use if damaged. If you are unable to connect to Bluetooth or to measure data, contact support from MedStorm.

Note.

For optimal functionality of the battery of the PainSensor, the PainSensor should not be stored fully charged or completely empty for a longer period of time. Then the battery lifetime may suffer. This is not of risk to the patient, however, optimal functioning of the battery is recommended.

5.1 Lifetime

The minimum lifetime of the system is 3 years on the condition that the instructions in this manual are followed. If the instructions of storing the device are followed, optimal functioning of the battery is sufficient for the lifetime of the device.

Preventive maintenance

The PainSensor does not need to be calibrated during the specified lifetime years, presuming the instructions in this manual are followed.

5.2 Support information

Email	support@med-storm.com
Tel	+47 907 88 976

5.3 Cleaning

Always disconnect the PainSensor from its power supply before cleaning. Ensure that the plug for protecting the charging USB port is closed before cleaning.

Note.

*Under no circumstances should the PainSensor and accessories be **immersed** in any liquid cleaning agent. Nor should it be exposed to steam or hot air sterilisation, or chemical sterilisation using ethylene oxide. Never use ether or petroleum-based solvents.*

For each time the device will no longer be used on the same patient, the electrodes used shall be removed, and the PainSensor and its accessories may be cleaned by wiping a clean cloth either with soap and water or with 70% isopropyl. For longer use on the same patient, the electrodes must be replaced after 48 hours.

5.4 Scrapping instructions

Dispose of this unit as normal electronic household waste per local regulations according to the WEEE Directive.

6. Appendix A - Set-up of the PSM connectivity box to connect with the PainSensor

The purpose of the PSM connectivity box is to receive data from the PainSensor and forward it to medical monitors, e.g. Philips IntelliVue Monitor. The data is transferred from the PainSensor to the PSM connectivity box using Bluetooth Low Energy, BLE.

Currently 2 versions are available. PSM connectivity box for Philips Monitors (#6015) need to be used to connect to Philips monitors and is shown in Figure 6-1. PSM connectivity box for Masimo Monitors (#6014) need to be used to connect to Masimo monitors. Read the instruction for use for these products in order to operate them.

Figure 6-1: Setup for connecting to Philips monitors with the PSM connectivity box

7. Appendix B - Connecting the PainSensor to a stand-alone PC with Bluetooth connection

The purpose of the PainSensor application is to receive data from the PainSensor and present it in the interface. The data is transferred from the PainSensor to the PainSensor application on a standalone PC / tablet using Bluetooth Low Energy, BLE.

Note.

Tablets from Med-Storm Innovation AS that are not medical approved shall not be placed in the patient zone (in the proximity of 1.5 meters from the patient).

7.1 The PainSensor Application

The interface of the PainSensor Application is shown in Figure 7-1.

Figure 7-1: The interface of the PainSensor Application

7.2 Connect to PainSensor

The PainSensor Application is connected to a PainSensor through bluetooth.

1. Switch on the PainSensor, and press “Connect Sensor” on the interface. The PS Application starts searching for PainSensors to connect with. When finding a PainSensor, the PS Application will connect to the sensor.

Note: If multiple sensors are found PC will connect to the sensor with the strongest signal.

2. When the PS Application is connected to the PainSensor and data is received, graphs showing the data will appear on the right side of the interface. Figure 7-2

Note.

Once the PainSensor is connected to a tablet delivered by MedStorm with PainSensor Software, this connection is remembered until you press “close” in the program. If you will continue to use the same sensor, simply press “disconnect sensor” after each patient, and when pressing “connect sensor” it connects to the same sensor.

Note.

When connecting for the first time, if multiple sensors are found, the software will connect to the sensor with the strongest signal.

Figure 7-2: The interface with example data

7.3 Description of the interface

In addition to the option of changing application modes, the different parts of the interface are described in this section.

7.3.1 Enter patient id

Entering the patient id gives you the option to save the measurement data to an ID. Do not use any personal data as a patient ID. You can either enter the patient ID before or during the measurement. If you forget to enter the patient ID before or during the measurement you will also have to option to enter this after clicking “Disconnect Sensor”.

7.3.2 Application modes

Press the field with “Anaesthesia” if you wish to change application mode. A drop-down menu will open, displaying the six different application modes that can be selected: Anaesthesia, PostOperative, ICU, Infants, WithDrawal and Neural Block, see image below. Select the preferred mode by pressing on it. Depending on which application mode is selected, the number of graphs on the right side of the application varies between 1-3. The different data that can be shown by the graphs are pain-nociceptive index, awakening index, nerve block index and withdrawal index.

- Anaesthesia: Pain-nociceptive index, awakening index and nerve block index
- PostOperative: Pain-nociceptive index
- ICU: Pain-nociceptive index and awakening index
- Infants: Pain-nociceptive index
- WithDrawal: Withdrawal/abstinence index
- Neural Block: Nerve block index

7.3.3 Raw data graph

This graph shows the raw data from the PainSensor. The y-axis shows the MicroSiemens level, whereas the x-axis shows the time. Note that the micro-Siemens baseline can differ between patients, which is normal.

7.3.4 Close, comment

The two buttons have different functions:

- “Close” should be pressed if you want to close the software and “forget” the PainSensor connection. Use this button if you wish to reconnect to the PainSensor another time or connect to a different PainSensor. If you want to connect to the same PainSensor next time you use the tablet, you can press disconnect instead, and the device will remember the PainSensor it last connected to.
- “Comment” can be used if you want to store some comment in the data. When you press “comment” you choose the moment you want to save the comment. That means, if you press comment at 09.03.00, the data will be stored at this time, independent of how long time you spend to write the comment.

7.3.5 Connect/disconnect to patient monitor and Connect/disconnect sensor

Connect/disconnect to patient monitor: Press this button in order to connect the tablet to a patient monitor and display the indexes on the patient monitors. Make sure you connect the tablet to a patient monitor beforehand. For example, to connect to a Philips patient monitor, insert the EC5 in the tablet that is connected to the EC10 with the Philips monitor. When pressing “connect monitor” ensure that the PainSensor is on and connected to the tablet beforehand.

Connect/disconnect sensor: Press “Connect sensor” to connect PainSensor to the PainSensor application. The measurement will automatically connect as long as the PainSensor is turned on and the LED is green. When you are done measuring, press “disconnect sensor”. A pop-up window will ask if you want to save the recording or not. You will also have the chance to enter patientID, if already added the patientID will appear. For more information about data storage, see section 5.

7.3.6 Zoom in / Zoom out

These two buttons give you the option of zooming either in or out on the trend graphs for PainIndex, Awakening Index and Nerve Block Index. The following options are available, 60 seconds, 60 minutes, 3 hours and 9 hours.

7.3.7 Pain Index

The PainIndex is displayed with a number and a graph. The graph is updated each second and based on the frequency of skin conductance peaks detected in the raw data. The shaded areas of the graph indicate the area where the PainIndex should be when a patient is receiving some stimuli. This is just an indication. However, it can be used to determine whether the patient would possibly need more or less analgesics. If the patient who is receiving analgesia is receiving a stimulus, and the graph is above the shaded area (pain-index of 4) it is an indication that this patient possibly needs more analgesia. If a patient who is receiving analgesia is receiving a stimulus and the graph is below 1, it is an indication that the patient has received too much analgesia. If a patient is not receiving any stimuli, an index of 0 does not suggest that the patient received too much. Once the patient has emerged from general sedation, it is recommended to switch to "post operative" setting in the application

Note.

Always use your own clinical judgment before changing levels of analgesia based on the information from the PainSensor.

7.3.8 Awakening index

The Awakening Index is displayed with a number and a graph. The graph is updated each second, and based on the size of skin conductance peaks detected in the raw data. If the graph goes above the shaded area (40), it is an indication that the patient is possibly about to wake up. Awakening index is a response that appears only if there is a strong nociceptive stimulus, it will not be triggered in case there is too low level of hypnotics, and the anti-nociceptive drug level is adequate.

7.3.9 Nerveblock index

The Nerveblock Index is displayed with a number and a graph. If the graph goes above the shaded area, it is an indication that the nerve is not blocked / no longer blocked.

7.4 Connect the tablet with a Philips patient monitor


1. Start the Philips patient monitor
2. Connect the Philips patient monitor by using an EC10 Intellibridge and EC5 (104) cable to VGA on the USB-C converter cable. The EC10 is without any light.
3. Launch the application on the tablet and connect to the sensor.
4. Wait until the sensor is connected and you can see the graphs.
5. Then take the USB-C cable into the tablet, and when the EC 10 is blinking green Click on the "Connect to monitor" button in the software application.
6. When the connection is established, the EC 10 turns green (if it turns blue it is something wrong, and restart the process).
7. After a couple of seconds, the indices should be displayed on the patient monitor.

7.5 Store and analyze data

7.5.1 Store files and locate them on the tablet

In order to store the data locally on the tablet:


- a) Enter patient id as explained in section 7.3.
- b) Press “disconnect sensor” after you are done with the measurement
- c) A pop-up will ask whether you want to “Save recording?”, enter patient data if not already added. Press “Yes” to save the data or press “No” to delete the data. It is also possible to save data without entering patient ID.



- d) Open “file explorer” in the menu found in the left corner.


- e) Press “Local Disk (C:)”


f) Press the folder “ProgramData”. If you cannot locate the folder “ProgramData”, press the bar on the top and type the name in.

g) Open the folder “PSS Application”

h) Your recordings are stored in this folder with the patientID, time and date included in the filename, there will be a raw data file stored as a .txt file and a .xlsx (excel) file. To analyze the data, follow the steps in the next section.

7.5.2 Analyse the files you stored

After locating the files on the tablet, insert a removable storage (either memory stick or and micro-SD card) into the tablet or computer, depending on the ports available. Copy the file you wish to analyse to the removable memory and open them on a computer that has Microsoft excel installed. Then open the files on this computer, and the data is shown. Here you will see the time, the indexes, and the comments you stored during the measurement.

8. Appendix C - Environmental and handling conditions

Measuring Unit		
Operating	Ambient temperature	+10°C – +40°C (50°F – 104°F)
	Ambient pressure	800hPa – 1060hPa (11.6 PSI – 15.4 PSI)
	Ambient humidity	30% - 75%
Transport	Ambient temperature	-10°C – +70°C
	Ambient pressure	500hPa – 1060hPa (7.3 PSI – 15.4 PSI)
	Ambient humidity	10% - 90%
Storage	Ambient temperature	+10°C – +30°C (50°F – 86°F)
	Ambient pressure	700hPa – 1060hPa (10.2 PSI – 15.4 PSI)
	Ambient humidity	30% - 75%
Transportation	It is possible to transport the system worldwide by air, road, ship and train.	
Altitude	The equipment cannot be used on altitudes higher than 2000 m above sea level.	
EMC/ESD	The PainSensor meets requirements in accordance with IEC 60601-1-2 Electromagnetic compatibility.	

9. Appendix D - Technical specifications

Measurements accuracy	Electronic noise which is below 0.01 microsiemens is filtered out, the electronic noise level is +/- 0.0025 µS (2SD). This applies for resistive measurements on 100 µS. PainIndex (peaks per sec): +/-0,13 AUC: +/- 87.11 microsiemens.sec
Measurement range	1-200 µS
Classification of medical device	Class II A
Maximum current definition	36 µA RMS The maximum value of current that can be supplied to a patient through the C electrode.
Power supply	The measuring unit operates on power from an external power supply of medical grade. Do not use any other than the provided power supply unit, unless it has been tested and verified by Med-Storm that it works together with the Measuring Unit Mains power input to measuring unit power supply is 90-264 VAC, 47-63Hz. Power consumption while loading is 2,5 W

Mechanical dimensions

Part	Weight [kg]	Dimensions [mm]
PainSensor	~0.08	85 x 70 x 22

List of cables and maximum lengths of cables

Cable	Maximum length [m]	Manufacturer	Model or part #
Mains cable, PainSensor	2	Med-Storm Innovation AS	with PainSensor (1100)

Hardware Specifications for the PainSensor Application (minimum requirements)

Operating System	Windows 10 or newer
Bluetooth	Bluetooth 5.x

Data security and network

A BLE connection based on BLE V5.0 is used. The transmission frequency is 2.402-2.480 GHz. Before the measured data of the patient is transmitted to a receiver, the receiver is verified by a verification algorithm. Only if the receiver passes the verification, the connection is maintained and the data is transmitted, otherwise the connection with the receiver is disconnected after verification. The transmitted data is encrypted by the Advanced Encryption Standard AES-128.

10. Appendix E - Regulations and Directives

The PainSensor meets applicable requirements of the following:

Reference	Title
MDR EU 2017/745	Medical Device Regulation (EU) 2017/745

GDPR (EU) 2016/679	The General Data Protection Regulation
WEEE Directive 2012/19/EU	Directive 2012/19/EU of the European Parliament and of the Council of 4 July 2012 on waste electrical and electronic equipment (WEEE)
RoHS Directive (EU) 2017/2102	Restriction of Hazardous Substances in electrical and electronic equipment. Directive 2017/2102/EU
RED Directive 2014/53/EU	The Radio Equipment Directive (RED), 2014/53/EU
CIR EU 2021/2226	Electronic instructions for use of medical devices

11. Appendix F - Electromagnetic compatibility

The PainSensor has been designed and tested to comply with the electromagnetic compatibility (EMC) limits for medical devices to the EN/IEC 60601-1-2 standard for EMC. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation.

The PainSensor does not need to be maintained or checked regarding EMC aspects during its service life.

WARNING:

Use of the PainSensor adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Portable HF communication devices (radios) (including their accessories such as antenna cables and external antennas) should not be used within 30 cm (or 12 inches) of the parts and leads of the PainSensor designated by the manufacturer. Non-compliance may lead to a reduction in the performance characteristics of the device.

CAUTION:

Use of accessories and cables other than those specified or provided by Med-Storm could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Medical electrical equipment requires special precautions regarding EMC and must be installed and operated according to these instructions. It is possible that high levels of radiated or conducted radiofrequency electromagnetic interference (EMI) from portable or mobile RF communications equipment or other strong or nearby radio-frequency sources, could result in performance disruption. Evidence of disruption may include erratic readings (e.g. strong fluctuations in the measuring signals), equipment ceasing to operate, or other incorrect functioning. If this occurs, survey the site to determine the source of disruption, and take the following actions to eliminate the source(s):

- Turn equipment in the vicinity off and on to isolate the source
- Relocate or re-orient the interfering source
- Increase distance between source and the PainSensor unit
- Educate clinical staff to recognize potential EMI related problems
- Restrict use of cell phones, etc. near the PainSensor unit
- Purchase medical devices that comply with IEC 60601-1-2

Guidance and manufacturer's declaration – electromagnetic emissions

The PainSensor is intended for use in the electromagnetic environment specified below. The customer or the user of the PainSensor should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The PainSensor uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The PainSensor is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Figure 11-1: Electromagnetic compatibility 201

The device complies with Part 15 of the FCC Rules

Guidance and manufacturer's declaration – electromagnetic immunity			
The PainSensor is intended for use in the electromagnetic environment specified below. The customer or the user of the PainSensor should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	+/- 8 kV contact +/- 15 kV air	+/- 8 kV contact +/- 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / Burst IEC 61000-4-4	+/- 2 kV for power supply lines	+/- 2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	+/- 0.5 kV, +/- 1 kV (Line-to-line)	+/- 0.5 kV, +/- 1 kV (Line-to-line)	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on	0 % UT for 0,5 cycle 0 % UT for 1 cycle	0% UT for 0,5 cycle 0 % UT for 1 cycle	Mains power quality should be that of a typical commercial or hospital environment. Use of the PainSensor during power mains interruptions is

power supply input lines IEC 61000-4-11	70 % UT for 25/30 cycles 0 % UT; 250/300 cycles	70 % UT for 25/30 cycles 0 % UT; 250/300 cycles	no problem, because the PainSensor is powered by a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE U _T is the a.c. mains voltage prior to application of the test level.			

Figure 11-2: Electromagnetic immunity 202

Guidance and manufacturer's declaration – electromagnetic immunity			
The PainSensor is intended for use in the electromagnetic environment specified below. The customer or the user of the PainSensor should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the PainSensor, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz 3 Vrms in the ISM bands	3 V	Recommended separation distance: $d = 1,2\sqrt{P}$ $d = 1,2\sqrt{P}$ 80 MHz to 800MHz $d = 2,3\sqrt{P}$ 800 MHz to 2,7 GHz
Radiated RF (IEC 60601-1-2:2014 4 th edition)	3 V/m 80MHz to 2700MHz	3 V/m	where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol.

NOTE 1 At 80MHz and 800MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflected from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the PainSensor is used exceeds the applicable RF compliance level above, the PainSensor should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the PainSensor.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Figure 11-3: Electromagnetic immunity 204

Guidance and manufacturer's declaration – immunity to high-frequency wireless communication equipment				
<u>Service</u>	<u>Test frequency in MHz</u>	<u>Frequency band in MHz</u>	<u>Modulation</u>	<u>Test level</u>
TETRA 400	385	380 – 390	18 Hz puls modulation	27 V/m
GMRS 460, FRS 460	450	430 – 470	1kHz sinus FM-modulated ± 5 kHz	28 V/m
LTE Band 13, 17	710, 745, 780	704 – 787	217 Hz puls modulation	9 V/m
GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	810, 870, 930	800 – 960	18 Hz puls modulation	28 V/m
GSM 1800, CDMA 1900, GSM 1900, DECT, LTE Band 1, 3, 4, 25; UMTS	1720, 1845, 1970	1700 – 1990	217 Hz puls modulation	28 V/m
Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	2450	2400 – 2570	217 Hz puls modulation	28 V/m
WLAN 802.11 a/n	5240, 5500, 5785	5100 – 5800	217 Hz puls modulation	9 V/m

Figure 11-4: immunity to high-frequency wireless communication equipment

12. Appendix G - Pre use checklist

Check	Signature
Verify that the system is not used on a patient with a skin condition which may affect skin conductance (e.g. injury of the skin).	
Verify that maximum one PainSensor is used on a patient with an electrically sensitive life support system (e.g. implantable pacemaker or defibrillator).	
Verify that the system is not going to be used when the patient has an injury affecting the sympathetic skin nerves .	
Verify that the system is not going to be used more than 48 hours in row, on the same patient, due to the electrodes that have to be changed.	
Verify that the electrodes are placed according to this manual. (Ch 4.2.1.1 for adults, Ch 4.2.1.2 for premature infants)	
Verify that the electrodes are of the correct type and approved by Med-Storm.	
Verify that if you temporarily disconnect any of the electrodes from the PainSensor, you see a magenta LED on the PainSensor.	
Verify that you have a secondary monitor to determine the patient's sensitivity to pain, such as e.g. blood pressure measurement.	
Verify that you are using MedStorms charger cable and charger plug unit to charge the device.	

For reasons of safety, the device may only be used if all the requirements above are satisfied.

PainSensor

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VERSIÓN 16

ESPAÑOL

DHF-00069

PainSensor REF 1002

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Manufacturer/Distributor:

Med-Storm Innovation AS
Gimle Terrasse 4
NO-0264 Oslo
Norway

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Référence / désignation	Description
MDR EU 2017/745	Medical Device Regulation (EU) 2017/745
GDPR (EU) 2016/679	The General Data Protection Regulation
WEEE Directive 2012/19/EU	Directive 2012/19/EU of the European Parliament and of the Council of 4 July 2012 on waste electrical and electronic equipment (WEEE)
RoHS Directive (EU) 2017/2102	Restriction of Hazardous Substances in electrical and electronic equipment. Directive 2017/2102/EU
CIR EU 2021/2226	Electronic instructions for use of medical devices

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Warnsymbol:	
	Caution is necessary when operating the PainSensor.
	Die Schachtel (und der Inhalt) sollten trocken gehalten werden
	Read the entire instructions for use before operating this PainSensor
	Class II electrical equipment for protection against electric shock (safety connection to electrical earth not required).
	Dispose of this unit as normal electronic household waste per local regulations according to the WEEE Directive
	TYPE BF APPLIED PART
	Serial number of device
	TYPE B APPLIED PART
	Product model identifier in catalogue
	Batch designation
	Date of manufacture
	Medical device
	Manufacturer
	Unique device Identifier
	Temperature limit
	Atmospheric pressure limitation
	Humidity limitation
	CE-Zeichen und Nummer der benannten Stelle
	QR-Code für den Online-Zugriff auf das Benutzerhandbuch

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Reference	Title
MDR EU 2017/745	Medical Device Regulation (EU) 2017/745
GDPR (EU) 2016/679	The General Data Protection Regulation
WEEE Directive 2012/19/EU	Directive 2012/19/EU of the European Parliament and of the Council of 4 July 2012 on waste electrical and electronic equipment (WEEE)
RoHS Directive (EU) 2017/2102	Restriction of Hazardous Substances in electrical and electronic equipment. Directive 2017/2102/EU
RED Directive 2014/53/EU	The Radio Equipment Directive (RED), 2014/53/EU
CIR EU 2021/2226	Electronic instructions for use of medical devices

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